KURARAY  
DENTAL

**Technical Manual**  
**Manuel technique**  
**Manual técnico**

***Epricord***<sup>TM</sup>

# **Contents**

***ENGLISH* ..... 1~31**

***FRANÇAIS* ..... 33~63**

***ESPAÑOL* ..... 65~95**

## Introduction

**Crown and bridge resins have been used extensively for restorations of crowns and bridges in routine practice because they are easy to handle and have relatively few contraindications. However, it is said that they are not completely satisfactory because there is room for improvement when comparing the esthetics and durability with porcelain.**

**EPRICORD™, which stands for Enhanced PProperties Integrated COmposite Resin Design, is our response to those shortcomings. Enhanced Properties refers to improved esthetics, handling and durability.**

**EPRICORD™ is our latest crown and bridge resin resulting from years of research and development into composite resins and adhesive materials. We believe we have achieved the highest level of handling ease and esthetic properties without compromising durability.**

## Contents

- Introduction
- Precautions
- Indications and contraindications
- Features
- Kit components
- Related products
- Basic laboratory work procedure
  1. Flow chart
  2. Light curing unit, polymerization time and depth of cure
  3. Color matching
  4. Fabrication of facing crown and bridge
  5. Fabrication of jacket crown
  6. Additional forming
  7. Repairing procedure
- Laboratory techniques
  1. Improving esthetics of the cervical area
  2. Producing natural transparency of the incisal edge
  3. Reproducing internal coloration like natural dentition  
(Cases where the thickness is extremely small)
  4. Color matching of colorless facing crown
  5. Combination with porcelain

## Enhanced PRoperties

**Integrated  
COmposite  
Resin  
Design**

### **Esthetic properties:**

Easy to achieve a high level of color matching

### **Easy handling:**

Superb fabrication efficiency

### **Durability:**

Hard to break and peel off

# Precautions

Cautions concerning the method of using the product

1. Use the light-curing unit and curing time specified in the "2. Light curing unit, polymerization time and depth of cure".
2. Do not mix different types of Body Resin together or Body Resin with other materials to prevent air entrapment in the resin or reduction in physical properties.
3. When loading the product, use caution not to expose it to strong light such as direct sunlight or use it immediately under laboratory light; otherwise the paste might harden.
4. The applicator cap should be replaced as soon as possible after resin has been dispensed from the bottle. This will prevent the evaporation of volatile contents, premature setting and contamination of the resin in the bottle.
5. If the area which needs polymerization is larger than the diameter of the tip of the visible light curing unit, it is advisable to divide the exposure session into a few installments.

## Caution to be taken when using the product

### 1. Warning

- 1) This product should only be used by dentists or dental technicians.
- 2) This product should not be used for purposes other than specified in Paragraph "Performance and indications".
- 3) When any of the related materials specified in these Instructions for Use is used with this product, read the Instruction for Use of the said material carefully.

Whether this product should be indicated or not must be determined by the dentist.

In particular, the product should not be used for cases with a covered occlusal surface of molars or close bites (eg. edge-to-edge bites) or bruxism.

### 2. Basic cautions to be observed

- 1) The product should not be used for patients with an allergy to methacrylate monomers.
- 2) If any hypersensitivity such as a rash or dermatitis occurs, discontinue the use of the product and consult a physician.
- 3) Avoid any direct contact with human body to prevent hypersensitivity. Wear gloves or take other appropriate measures when using the product. In particular, avoid touching unpolymerized resin.
- 4) Use caution to prevent the product from coming in contact with the oral soft tissue or skin or getting into the eye. If the product comes in contact with the oral soft tissue or skin, wipe it off with a cotton pledget or a piece of gauze moistened with alcohol and immediately wash with copious amounts of water. If the product gets in the eye, immediately wash the eye with copious amounts of water and consult a physician.
- 5) Use caution to prevent someone from swallowing the product accidentally. Avoid using the same preparation for different patients to prevent cross infection.
- 6) Use safety equipment such as a local vacuum unit, protective shield, and dustproof mask approved by an appropriate organization, to prevent the operator from inhaling dust when shaping, finishing or glazing the product.

## Storage and Expiry Date

### Storage

- The product should not be stored in direct sunlight. It should be stored at room temperature (4°C -25°C /39°F -77°F).
- Opaque Primer is flammable. Do not use or store it near open flame.

### Expiry date

The product should be used by the date indicated on the package. \*  
The expiry date is determined by self-certified (using data collected by Kuraray Medical Inc.)

\* Exp. 2002-05 means that the expiry date of the product is May 2002.

# Indications and contraindications

## **Indications**

This product is designed to be used for the following indications:

- Crown restorations using facing crowns, facing bridges or jacket crowns
- Repairs of the crown and bridge resin or porcelain of facing crowns and jacket crowns

## **Contraindications**

- Patients with an allergy to methacrylate monomer
- Covering the occlusal surface of molars
- Close bites (eg. edge-to-edge bites) or bruxism

# Features

## 1. Improved color matching

Epicord™ reproduces accurate color matching of the patient's original tooth color in accordance with the Vita shade guide (classical). With basic colors reviewed according to the Vita shade guide (classical), it reproduces vivid texture like natural tooth structure by the double-layer method.

### 1) Excellent color matching with surrounding natural teeth

There is less chance of color having a yellowish or whitish tint in the oral cavity, a phenomenon often observed with conventional crown and bridge resins when they are loaded. With the double-layer method, superb color matching with the surrounding natural teeth can be achieved.

21] : Epicord™

### 2) Improved color matching ability on thin layers

It is said that color matching is difficult with facing crowns because the application layer of crown and bridge resins is thinner than porcelain.

Epicord™ offers exceptional color matching thanks to its strong color development by applying Enamel in double layers.

Comparison with Vita shade guide(1 : Epicord™)

## 2. A high level of color matching attained by simple operation

A high level of color matching can be achieved by simple operation using just two Opaque resins, Cervical and Incisal.

### 1) Cervical Opaque

Cervical Opaque is used for cases where the Opaque color reflects strongly or for areas close to the cervical area to offer more exact color matching. In addition, it may be used for adding a reddish tint to the A type color of the Vita shade guide (classical).

Applied(1|1)

Not applied(1|1)

### 2) Incisal Opaque

Incisal Opaque offers transparency similar to that of natural teeth for cases where a metal backing is required at the incisal edge.

### 3. Increased number of shades available for a variety of color matching

The number of color variations has been increased to make it possible to offer a wide range of color matching without requiring high level of skill, as often needed when using porcelain.

#### 1) NP shades

It is said that almost all natural tooth colors come under Type A of the Vita Shade Guide. However, Vita shade A series do not necessarily correspond to all the shades of natural teeth. Epicord™ has new additions its component kit, NP1.5 and 2.5 shades (pinkish orange) which are not covered in the Vita Shade Guide (Classical), to make it possible to be indicated for a wider range of cases than before.

11 : Epicord™

#### 2) NW shades

Keeping pace with the recent trend of teeth whitening treatment, NW0 and NW0.5 have been added to respond to the demand of patients for whiter teeth as well as to maintain the color after whitening.

111 : Epicord™

#### 3) Opacious Dentins

Opacious Dentins are used when the build up is ground thin or when a colorless facing crown or the entire or part of a pontic is repaired using resin.


#### 4) Effects

Five types of Effect colors are available for more authentic color matching at incisal edges.

#### 5) Other colors


There are two Gingival colors and two Opaque Modifier colors, W (white) and P (pink).

## 4. Economical product characteristics

### 1) Improved shaping ability

The paste is not too hard yet maintains its form due to proper viscosity, making it easy to contour prominences on tooth surfaces.

### 2) Polishability

Polishability has been improved while maintaining excellent wear resistance and surface smoothness.

### 3) Short curing time

The preliminary polymerization time for Body Resin has been reduced to improve work efficiency substantially.

### 4) Use of half-size syringes

The quantity of shades not used frequently has been reduced in half to alleviate purchasing and inventory costs.


## 5. High strength and adhesiveness to metal

Epicord™ offers high physical properties in addition to esthetic properties and easy handling. Needless to say about the strength and durability of Body Resin, it combines with a metal-adhesive primer containing phosphate monomer (MDP) to provide a strong bond for the resin to a metal frame, producing an exceptionally esthetic restorations.


1: Epicord™ 1: Porcelain


### Compressive strength

### Flexural strength


### Surface hardness


### Metal adhesiveness

Adhesive strength of Opaque Resin to gold-silver-palladium alloy (after 30-day immersion in water)

### Wear resistance

(The amount of wear of tooth brush)

## Kit components

### Opaque Resin

This is a single-paste light-cured composite resin consisting of a matrix of multifunctional methacrylate monomer, designed for masking out metal color.

Opaque paste has optimal viscosity that allows easy application with a brush. There is less development of unpolymerized resin on the surface after light curing. It is chemically cured on contact with Opaque Primer to provide a strong bond for the Opaque layer to a metal frame.

### Body Resin

Body Resin is a light-cured, highly filled hybrid composite resin that contains ternary fillers (particulate inorganic filler, inorganic microfiller, and organic composite filler) in a matrix of multifunctional methacrylate monomer.

### Opaque Primer

Opaque Primer contains metal-adhesive phosphate monomer (MDP) to provide strong adhesion of the Opaque resin to the metal frame.

### Repair Liquid

Repair Liquid is a low-viscosity, light-cured resin used when shaping the paste or when loading additional loading of paste.

## System components of Epricord™

Component	Standard set	A3 set	Components separately sold (The following components are not included in each set and separately sold)
<b>1. Opaque Resin</b>			
① Cervical Opaque	CO1, CO2, CO3: one for each (2.0ml)	CO1:one (2.0ml)	CO4, CO5: (2.0ml)
② Body Opaque	OA1, OA2, OA3, OA3.5, OA4, OB2, OB3, OC3: one for each (2.0ml)	OA3:one (2.0ml)	OB1, OB4, OC1, OC2, OC4, OD2, OD3, OD4 (2.0ml) ONW0,ONW0.5 ONP1.5,ONP2.5 (2.0ml)
③ Incisal Opaque	IO1, IO2:one for each (2.0ml)	IO1:one (2.0ml)	—
④ Opaque Modifier	—	—	W,P (1.0ml)
<b>2. Body Resin</b>			
① Cervical	CE1, CE2, CE3, CE4, CE5: one for each (3.8g/2.6ml)	CE1:one (3.8g/2.6ml)	CE6, CE7, CE8 (3.8g/2.6ml)
② Opacious Dentin	ODA1, ODA2, ODA3, ODA3.5, ODA4: one for each (3.8g/2.6ml)	ODA3:one (3.8g/2.6ml)	—
③ Dentin	DA1, DA2, DA3, DA3.5, DA4, DB2, DB3, DC3: one for each (3.8g/2.6ml)	DA3:one (3.8g/2.6ml)	DB1, DB4, DC1, DC2, DC4 DD2, DD3, DD4 (3.8g/2.6ml) DNW0,DNW0.5, DNP1.5,DNP2.5 (3.8g/2.6ml)
④ Enamel	E1, E2, E3:one for each (3.8g/2.6ml)	E1:one (3.8g/2.6ml)	E0 (3.8g/2.6ml) E4 (3.8g/2.6ml)
⑤ Transparent	T:one (3.8g/2.6ml)	T:one (3.8g/2.6ml)	TLV (3.8g/2.6ml)
⑥ Effect	—	—	CT, AM, WE, TB, CE (3.8g/2.6ml)
⑦ Gingival	—	—	P1, P2 (3.8g/2.6ml)
⑧ Base Material	one (3.8g/2.8ml)	—	—
<b>3. Opaque Primer</b>	one (9ml)	one (9ml)	—
<b>4. Repair Liquid</b>	one (6ml)	one (6ml)	—
<b>5. Accessories</b>			
① Flat brush (#933)	2pcs.	1pc.	—
② Small brush (#901)	2pcs.	1pc.	—
③ Small brush (#932)	2pcs.	1pc.	—
④ Mixing plate (#954)	8pcs.	1pc.	—

All the components of each set given above are separately available as refills.

## Related products

Metal-adhesive primer

### Alloy Primer

Alloy Primer is a metal-adhesive primer composed of cationic adhesive monomer and phosphate adhesive monomer (MDP). The application of this primer to precious metal alloys offers the same level of adhesion enhancement effect as tin-plating.

Separating agent for jacket crowns

### Jacket Separate Kit

A plaster separating agent for jacket crowns. It consists of Jacket Spacer to provide a cement space and Margin Sep to separate plaster from the jacket crown. Use of this kit allows the jacket crown to easily release from the working model.

Light-cured staining material

### ChromaZone™ Color Stain

A light-cured, low-viscosity composite resin that is used for the characterization and color correction of the crown.

Dental composite resin bonding agent

### Clearfil™ Repair

The bonding material that is used for intra-oral repair with Epricord™.


# Basic laboratory work procedure

## 1. Flow chart

The following explanation is given, taking the fabrication of a facing crown and a jacket crown as examples.

## 2. Light-curing units and polymerization time

(sec.)

	Opaque Resin	Body Resin	
		Preliminary Polymerization	Final Polymerization
UniXS	180	30	240

### 3. Color matching

#### List of shades

	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4	NW 0	NW 0.5	NP 1.5	NP 2.5			
<b>Body Resin</b>																							
Cervical	—	CE1		CE2		—	CE3		CE4		—	CE5		CE6		CE7		CE8		—	—	—	—
Opacious Dentin	ODA1	ODA2	ODA3	ODA3.5	ODA4	—	—	—	—	—	—	—	—	—	—	—	—	—	—	—	—	—	
Dentin	DA1	DA2	DA3	DA3.5	DA4	DB1	DB2	DB3	DB4	DC1	DC2	DC3	DC4	DD2	DD3	DD4	DNW 0	DNW 0.5	DNP 1.5	DNP 2.5			
Enamel	E1		E2		E4	E1	E3		E1	E3		E1	E3		E1	E3		E0	E1				
Transparent	Y (Transparent), TLV (Transparent Low Value)																						
Effect	TB (Transparent Blue), WE (White Enamel), CT (Cervical Transparent), CE (Creamy Enamel), AM (Amber)																						
Gingival	P1, P2																						
Base Material	BASE MATERIAL																						
<b>Opaque Resin</b>																							
Body Opaque	OA1	OA2	OA3	OA3.5	OA4	OB1	OB2	OB3	OB4	OC1	OC2	OC3	OC4	OD2	OD3	OD4	ONW 0	ONW 0.5	ONP 1.5	ONP 2.5			
Cervical Opaque	—	CO1				—	CO2				—	CO3		CO4	CO5	—	—	—	—				
Incisal Opaque	IO1, IO2																						
Opaque Modifier	W (White), P (Pink)																						

#### Applications of Body Resin

- Cervical : Paste used for producing color at cervical areas.
- Dentin : Paste used for producing dentin color.
- Enamel : Paste used for producing enamel color.
- Transparent : Paste used for producing transparency at incisal edges.

T(Transparent) : Producing usual transparency.  
 TLV(Transparent Low Value) : Transparent with low luminosity.

- Effect : Characterizing colors.

CT (Cervical Transparent) : Transparent used for producing transparency at cervical areas.  
 AM (Amber) : Producing amber enamel color.  
 WE (White Enamel) : Producing cloudy enamel color seen at the angular areas of anterior teeth and occlusal areas of premolars.  
 CE (Creamy Enamel) : Producing enamel color of low transparency.  
 TB (Transparent Blue) : Producing bluish enamel color.

- Opacious Dentin : Paste used for cases when the build up is ground thin or when a colorless facing crown or the entire or part of a pontic is repaired using resin.
- Gingival : Paste used for producing gingival color.

P1 ...Gingival color paste of high luminosity.  
 P2 ...Gingival color paste of low luminosity.


- Base Material : Paste used for filling up pontic holes.

#### Applications of Opaque Resin

- Body Opaque : Masking out metal color and producing basic colors of the crown.
- Cervical Opaque : Used when the thickness at cervical areas is small.
- Incisal Opaque : Producing transparency at incisal edges (two colors are available; IO1, IO2).
- Opaque Modifier : Adjusting Opaque Resin colors or using independently/ Two colors are available; P (pink) and W (white).

## Color matching procedure


### · Basic method (double-layer method)

· **Dentin** : Form the Dentin paste on the crown to shape the area from the cervical edge to 2/3 of the crown length into the same as the final contour. Make the paste thickness at the remaining 1/3 area gradually thinly toward the incisal edge.

· **Enamel** : Form the Enamel paste thinly on the crown from the incisal edge to 1/3 of the crown length. Form the Enamel paste on the proximal areas from incisal edge to 1/2 of the crown length to producing internal coloration.

### · Quadruple-layer method

· **Cervical** : Form the Cervical paste thinly on the crown from the cervical area to 1/4 to 1/5 of the crown length, making the thickness gradually thinner from the margins to the cervical area.

· **Dentin** : Form the Dentin paste on the crown to shape the area from the cervical edge to 1/2 to 1/3 of the crown length into the same as the final contour. Make the paste thickness at the remaining 1/2 to 1/3 area gradually thinly toward the incisal edge.

· **Enamel** : Form the Enamel paste thinly on the crown from the incisal edge to 1/5 to 1/6 of the crown length.


· **Transparent** : Form the Transparent paste thinly on the crown from the incisal edge to 1/2 to 1/3 of the crown length.

## 4. Fabrication of facing crowns and facing bridges

The product should be polymerized referring to 2. Light curing units and polymerization time on page 14.

### 1. Fabricating a working model

Fabricate a working model in the usual manner.

### 2. Fabricating a metal frame

#### 1) Waxing up

Apply a wax separating agent on the working model and form the shape of a crown on the model.

#### 2) Cutback

Cut back the area to be faced. Epicord™ has excellent strength, wear resistance, and metal adhesiveness so that the incisal edges and proximal areas are restored with the resin without the coverage of metal.

#### 3) Applying retention beads

Apply an adhesive on the facing surface and sprinkle retention beads to complete the wax pattern.

Use spherical retention beads with a diameter of 100µm-200µm. Improper use of retention beads causes exfoliation of the metal frame.

#### 4) Completing metal frame

Cast a metal frame in the usual manner. Check the casting for adaptation and finish it.

### 3. Sandblasting facing surface

Sandblast the facing surface of the metal frame using 50µm alumina powder. Clean the surface with ultrasonic cleaner or using a steam cleaner before drying.

### 4. Treating the facing surface of precious metal alloy

To enhance adhesion to metal, apply a metal-adhesive primer (eg. Alloy Primer) on the facing surface of the metal frame with the attached small brush and let it dry at room temperature. Or use a dental tin plating unit to tin plate the facing surface of the metal frame.

### 5. Application of Opaque Resin

It is indispensable to apply Opaque Primer before using Opaque Resin in order to cure Opaque Resin completely as well as to provide strong bond for Opaque Resin to metal.

#### 1) Applying and drying Opaque Primer

Apply Opaque Primer on the facing surface of the metal frame with the attached small brush and evaporate volatile content by blowing air lightly or leave it for 30 minutes.

#### 2) Applying and curing the 1st layer of Opaque Resin


Apply the Opaque Resin paste (eg. Body Opaque or Cervical Opaque) of the shade you selected, on the facing surface of the metal frame thinly with the attached flat brush. After application, cure the resin for the specified time.

Multifunctional methacrylate, the main component of Opaque Resin and Body Resin, has anaerobic setting characteristics, which keep the resin surface from hardening so long as it is exposed to the air. Do not remove this unpolymerized layer because it serves as a bond to the resin to be loaded next.

### 3) Filling up pontic hole with Base Material

When fabricating a facing bridge, fill Base Material to the pontic area for leveling the facing space. Cure Base Material for the specified time.

### 4) Applying and curing the 2nd layer of Opaque Resin

Apply Opaque Primer on the cured 1st layer of Opaque Resin and Base Material and dry it. Apply another layer of Opaque Resin on it and cure for the specified time. Repeat this procedure until the metal color is masked out. Instead of Opaque Primer and Opaque Resin of Epricord™, you can use Cesead™ II Opaque Resin and Opaque Primer. (However, the colors of Epricord™ may not be reproduced properly if Cesead™ II Opaque Resin is used.)

### 6. Forming and polymerization of Body Resin

#### 1) Forming and preliminary polymerization of Cervical

Form the Cervical paste thinly in the area from the cervical area to 1/4 to 1/5 of the crown length, making the thickness gradually thinner from the margins to the cervical area. Cure the paste preliminarily.

#### 2) Forming and preliminary polymerization of Dentin

To produce the basic color of the crown, form a generous amount of the Dentin paste and cure preliminarily.

#### 3) Forming and preliminary polymerization of Enamel

Apply a thin coat of the Enamel paste in the area from the incisal edge to 1/5 to 1/6 of the crown length and cure preliminarily.

### 4) Forming Transparent

Apply a thin coat of the Transparent paste in the area from the incisal edge to 1/2 to 1/3 of the crown length.

### 5) Final polymerization

After forming the Transparent paste in the contact area, perform final polymerization.

### 7. Shaping

Use a carbide bur, carborundum point, or silicon point.


### 8. Finishing

Using a silicon point, shape into a detail anatomy and smoothen the boundary between the resin and the metal while removing the scratch impressions due to the carborundum from the resin surface.

### 9. Polishing

Polish the resin using a brush or cloth buff with fine-particle polishing agents. Instead of polishing process, a dental resin surface treatment agent may be used. (However, in this case, the luster will disappear in the long run.)

### 10.

## 5. Fabrication of jacket crowns

The product should be polymerized referring to 2. Light curing units and polymerization time on page 14.

### 1. Fabricating a working model

Divide and trim the die and complete a working model. Use caution not to give too deep an undercut below the margins when trimming; otherwise the plaster model might break when releasing the crown from it. If you cut too deep an undercut, block out with wax.

### 2. Application of resin separating agent (eg. Jacket Separate Kit)

#### 1) Applying and drying Jacket Spacer

To make it easy to release the jacket crown from the die, apply Jacket Spacer evenly on all the area except the margins of the abutment. Dry the applied Jacket Spacer by leaving it or by blowing air gently. Apply Jacket Spacer repeatedly as necessary.

#### 2) Applying and drying Margin Sep

Apply Margin Sep on the margins of the abutment and dry by leaving it or by blowing air gently.

### 3. Application and polymerization of Opaque Resin

Apply Opaque Primer on the abutment and dry it. Apply Opaque Resin and cure for the specified time. If there is no need to mask out the abutment color (such as when abutment is tooth structure or composite resin), you can skip this step.


### 4. Forming and preliminary polymerization of Dentin

Refer to 6. Forming and polymerization of Body Resin of 4. Fabrication of facing crowns and facing bridges on page 17.

### 5. Forming and preliminary polymerization of Enamel and Transparent Refer to 6.

Forming and polymerization of Body Resin of 4. Fabrication of facing crowns and facing bridges on page 17.

### 6. Releasing of jacket crown from the die and checking for adaptation

Release the jacket crown from the die, remove resin separating agent adhering to the abutment or the inner surface of the jacket crown. Check for adaptation.

### 7. Shaping and polishing

Shape and polish the jacket crown to complete it.

## 6. Additional forming

If an additional forming or reshaping is needed after shaping, ground the area where you want to add paste to create a rough surface, apply Repair Liquid, and load paste before reshaping.

The polymerization should be performed, referring to 2. Light curing units and polymerization time on page 14.

### 1. Expose a fresh surface where an additional forming is to be made

Using a carborundum point, cut the area where you want to add paste to expose a fresh surface. Remove cuttings by blowing air.

### 2. Apply Repair Liquid

Apply a thin coat of Repair Liquid on the fresh surface. Instead of Repair Liquid, you can use Estenia™ Modeling Liquid.

### 3. Form Body Resin additionally

Form the Dentin paste on the dentin area and the Enamel paste on the enamel area. Cure the resins.

### 4. Shaping, polishing and finishing

Shape a rough contour using a carbide bur, carborundum point, or silicon point. Then shape to a detail anatomy while removing the scratch impressions due to the carborundum from the resin surface. Finally, finish the restoration by polishing.

## 7. Repairing procedure


When the dental visible light curing unit is used for intra-oral repairs, it is advisable to divide the exposure session into a few installments according to the area to be polymerized.

The product should be polymerized referring to 2. Light curing units and polymerization time on page 14.

### Cut fractured surface

Cut the fractured surface by one layer, using a diamond point.  
Give bevels on the incisal edges.  
(If the metal is exposed, roughen its surface.)

### Acid etch the fractured surface

Apply a phosphoric acid etching agent (eg. K-Etchant Gel) on the adherent surface of the facing or jacket crown you are going to repair, and leave it for 5 seconds before washing and drying.

When the metal is exposed to view

### Apply a metal-adhesive primer

Apply a metal-adhesive primer (eg. Alloy Primer) on the metal surface with a small brush and let it dry at room temperature.

### Apply a bonding agent

Treat the adherent surface to be repaired using a dental bonding agent for porcelain adhesion (eg. Clearfil™ Repair). For how to use, see the Instructions for use of the said bonding agent.

When the metal surface is exposed to view

### Apply and cure Opaque Resin

If the metal surface is exposed to view, apply a thin coat of the resin and cure, and then repeat the same procedure again to mask out the metal color.

### Form and cure Body Resin

Form the Body Resin pastes (Dentin, Enamel, Transparent, Gingival) and light cure.

### Shaping, finishing and polishing

Shape the restoration using a fine-particle diamond point or carbide bur. After shaping, polish using a white point and silicone point in turns.


# Laboratory techniques

The product should be polymerized referring to 2. Light curing unit and polymerization time on page 14.

## 1 Improving esthetics at cervical areas

If the Opaque color reflects so strongly at the cervical area that the color is whitish, use Cervical Opaque to improve color reproduction.

Cervical Opaque is applied thinly over Body Opaque. (If Cervical Opaque is directly applied on the cervical area in stead of Body Opaque, the color may become too dark.)

### 1. Apply and cure Body Opaque

Apply Body Opaque of a specified shade to produce a primary coat.

### 2. Apply and cure Cervical Opaque

Apply a thin coat of Cervical Opaque 2 to 3 mm in width from the margins to the cervical area, to the extent that you can see the color of Body Opaque through the Cervical Opaque.

### 3. Form and cure Dentin

#### Form Dentin

(Remember that if Cervical Opaque is used, there is no need to use Cervical Paste.)


### 4. Form and cure Enamel

### 5. Shape and finish


Comparison of cases when Cervical Opaque is used and when Cervical Opaque is not used

Repair examples of cases when Cervical Opaque is used and when Cervical Opaque is not used.

### Color matching for other shades

The procedure for using Cervical Opaque is the same as shown above.


After application of Opaque

Completed

## 2 Adding reddish color

To reproduce more reddish color than Vita Classical shades, it is advisable to apply a thin coat of Cervical Opaque in the area from the margins to the middle of the crown to the extent that the color of the prime coat can be seen through. Color adjustment using Opaque makes it possible to reproduce the color you want with a thin application.


### 1. Application and polymerization of Body Opaque (Completion of prime coat)

Apply Body Opaque of a specified color. The prime coat is completed using Body Opaque alone.

### 2. Application of Cervical Opaque

Apply a thin coat of Cervical Opaque in the area from the margins to the middle of the crown to the extent that the color of the Body Opaque can be seen through.


### 3. Forming Body Resin

Form Dentin and Enamel.

### 4. Completed

#### Repair examples


Colors obtained when Opaque is thinly applied

Color samples

### 3 Producing natural transparency at incisal edges

When a metal backing is required at incisal edges, transparency like natural dentition can be reproduced by using Incisal Opaque.

Build up diagram

#### 1. Apply and cure Body Opaque

Apply and cure Body Opaque in the usual manner.

#### 2. Apply and cure Incisal Opaque

Apply Incisal Opaque in such a way that prominences may be reproduced on incisal edges. A real transparency can be reproduced.

#### 3. Form Dentin

Apply a thin coat of Dentin at incisal edges.

#### 4. Form and cure Enamel

The photo shows Enamel formed only on the upper right central incisor.


#### 5. Completed

#### Repair example ( 111 )

## 4 Restoration using NP shades

Use of NP shades has made it possible to cope with pinkish orange tooth colors that could not be covered by the conventional Vita Shade Guide (Classical), thus offering improved color matching.

**1. Polymerize Opaque Resin and Cervical**

**2. Form Dentin**

**3. Apply Stain**

**4. Form Enamel**


**5. Completed**

**Repair example** (1. Upper right central incisor: Comparison of Epricord™ and NP shade guide)

## 5 Restoration using NW shades

Cosmetic teeth whitening such as bleaching has become popular. Epricord™ has new additions, NW shades, to its component kit to respond to the demand of patients for whiter teeth as well as to maintain the color after whitening. NW shades (NW0 and NW0.5) offer higher brightness than Vita Shade Guide (Classical).

**1. Form and cure Opacious Dentin**

**2. Forming Dentin**  
Form DA1 in the area from the cervical region to the cutting edge.

**3. Forming Dentin**  
Form DNW0.5 in the area from the cutting edge to the middle of the crown.

**4. Form Enamel**

**5. Completed**


**Repair example (Upper right and left central incisors) ( 11 )**

## 6 Producing internal coloration like natural dentition

(Cases where the thickness is extremely small.)

When the thickness is extremely small and the Opaque color appears strongly, use Opacious Dentin to produce internal coloration.

1. Working model

2. Apply and cure Opaque


3. Form Opacious Dentin

Make the thickness of the build up thin (about 0.15 mm) to create a scaled-down contour of the crown.

4. Form Dentin

5. Form Enamel


6. Completed

Comparison with Vita shade guide (classical)

## 7 Color matching of colorless facing crown

Esthetics as cervical areas as well as at incisal edges can be improved by fabricating a colorless facing crown.

### 1. Preparation form of abutment teeth

Prepare a shoulder at cervical area at the labial side to have a sufficient thickness.

### 2. Apply and cure Opaque Resin

Apply Opaque Resin at the facing area in such a way that the plaster model of the cervical area is covered slightly.

### 3. Form the cervical area

Form Opacious Dentin or Cervical on the crown from the cervical area to 1/4 to 1/5 of the crown length, making the thickness gradually thinner from the margins to the cervical area.

### 4. Form Dentin

Apply Chroma Zone™ Color Stain on the incisal edges.

### 5. Form Enamel and Transparent.

### 6. Shape and finish.

Repair example (upper left central and lateral incisors) ( 21 )

## 8 Combination with porcelain

The colors of Epricord™ are developed based on Vita shade guide (classical).

However, the A type colors are close to those of Noritake Super Porcelain AAA.

Therefore, it is easy to perform color matching when using both or adding to porcelain already installed.


**1. Form Cervical**

**2. Form Dentin**

**3. Completed**

### Repair example

(Upper left central incisor: porcelain, upper right central incisor: Epricord™)

## 9 Forming paste, taking advantage of physical properties of the material.

The Body Resin of Epicord™ is not too hard yet maintains its form making it easy to contour prominences on tooth surfaces.

In addition, it has less development of unpolymerized layer\* so that there is not need to build up a greater amount of paste to allow for the development of unpolymerized resin, offering efficient shaping.

Here are some points to be kept in mind when building up the Body Resin pastes.

\* The unpolymerized layer means a soft layer on the resin surface remained after the resin is cured.

### Form Dentin


#### 1.Placing the paste and holding it under pressure

Put the paste on the frame and disperse it lightly while forcing out air from inside the resin.

#### 2.Checking transitional area from metal

Finish the restoration smoothly using an instrument or a small brush to seal the transitional area between metal and the paste securely.

#### 3.Shaping

Reproduce prominences using an instrument with Repair Liquid applied to its tip or using a small brush.

#### 4.Checking margins

Check the sealing effect at the cervical region before polymerization.

### Form Enamel

#### 1.Placing the paste and hold it under pressure

There is no need to hold the paste under strong pressure when forming it.


#### (Use of Repair Liquid)

It is good practice to apply a thin coat of Repair Liquid to the surface of the cured resin to provide better wetting with the paste you are going to form.

#### Use of Repair Liquid is effective in the following case:

- Complicated prominences were given.
- Enamel paste was dispersed thinly.
- Resins were formed into multiple layers. (when Transparent and Effect were used.)

#### 2.Finishing

Finally, smoothen the transitional area between the cured Dentin and Enamel pastes using a small brush. This makes it efficient to shape and polish the restoration.

#### 3.Completion of forming

It is possible to form into almost the same size as the final restoration because there is less development of unpolymerized resin.

[Note]

CLEARFIL is a trademark of KURARAY CO., LTD.

[Warranty]

KURARAY MEDICAL INC. will replace any product that is proved to be defective.

KURARAY MEDICAL INC.

does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

Manufactured by  
**KURARAY MEDICAL INC.**  
1621 Sakazu, Kurashiki, Okayama 710-0801, Japan

European Union Representative  
**KURARAY EUROPE GmbH**  
Building F821, Hoechst Industrial Park  
65926 Frankfurt am Main, Germany  
Tel:+49 (0)69 305 35 840 Fax:+49 (0)69 305 35 640  
I2111EU1-02 Printed in Japan 10/2006

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The logo for Epricord, featuring the brand name in a bold, italicized, white sans-serif font. A red swoosh underline is positioned above the letters 'i', 'r', and 'c'. A small 'TM' trademark symbol is located at the bottom right of the word.

**Epricord**<sup>TM</sup>

## Enhanced PRoperties

**Integrated  
COmposite  
Resin  
Design**

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CLEARFIL es una marca registrada de KURARAY CO., LTD.