POLYMER-BASED CROWN AND BRIDGE MATERIAL

ENGLISH INSTRUCTIONS FOR USE

I. INTRODUCTION
ESTENIA™ C&B is a polymer-based crown and bridge material for creation of facing cast crowns, facing cast bridges, jacket crowns, inlays, onlays and bridges with frameworks.

II. INDICATIONS
ESTENIA™ C&B is indicated for the following applications for the restoring crowns and defects.

- 1) Facing cast crowns and facing cast bridges
- 2) Jacket crowns
- 3) Inlays and onlays
- 4) Bridges with frameworks

III. CONTRAINDICATION
This product should not be used for patients with a history of hypersensitivity, such as rash or dermatitis, when exposed to methacrylate monomer.

IV. PRECAUTIONS
Please review this Instructions for Use carefully before using the product. After reviewing the Instructions for Use, keep them in a convenient place for quick reference if necessary. If this Instructions for use are lost, contact the distributor or supplier shown on the package.
This product should only be used by dentists or dental technicians.
To prevent the occurrence of trouble with operation or storing, the product should be used according to the instructions specified in this Instructions for Use. Do not use this product for contraindications given in this Instructions for Use.

Safety precautions
The use of this product should be determined by a dentist individually for each patient.
This product should NOT be used when the following conditions are present:
- Malocclusion, clenching or bruxism
- Closed bite or cases where there is no occlusal support for biting except the tooth being treated (if used in these situations, the occlusal contact area should be covered with metal.)
- The need to cover the occlusal surfaces of molars (The product can be used if there is some occlusal support other than the tooth under treatment.)
- Clasp teeth for partial dentures (The tooth area which comes in contact with the clasp should be covered with metal.)
- Patients requiring bridges replacing more than two missing teeth, or with an edentulous span of more than 15 mm.
If any hypersensitivity, such as rash or dermatitis occurs, discontinue the use of the product and get medical attention.
Avoid any direct contact of material components with the human body to prevent hypersensitivity.
Wear gloves or take other appropriate measures when using the product. In particular, avoid touching unhardened resin.
Use caution to prevent the product from coming in contact with the oral soft tissues or the skin or getting into the eyes. If the product comes in contact with oral tissue or the skin, wipe it off with a cotton pledget or a piece of gauze moistened with alcohol and then wash immediately with copious amounts of water. If the products get into an eye, immediately wash the eye with copious amounts of water and consult an ophthalmologist.
Use caution to prevent anyone from accidentally swallowing the product.

Use safety equipment, such as a local vacuum unit, protective shields, and an appropriately approved dustproof mask, to prevent the operator from inhaling dust when shaping, finishing or polishing the product.
When any of the related materials covered in this manual are used with this product, read the Instructions for Use for that material carefully.
Handling and manipulation precautions
1) Avoid cracked or broken crowns or inlays
- Assume that the crown or inlay is sufficiently thick, using the methods for preparing abutments and cavities described in this manual.
- Assume correct curing by using the light curing units and heat curing units as specified in this manual. Observe the specified curing conditions.
- Review the part of additional fabrication procedures in this manual for additional fabrication.
- Do not mix different types of Body resin together, or mix Body resin with other materials, in order to prevent impairing the physical properties of the materials.
2) Avoid broken or cracked bridges with frameworks
- Assume that the bridge is sufficiently thick, using the methods for preparing abutments and cavities described in this manual. Take particular care not to use this product in cases where sufficient strength cannot be obtained.
- Assume correct curing by using the light curing units and heat curing units as specified in this manual. Observe the specified curing conditions.
- To prevent the impairment of product physical properties, do not use framework or crown materials supplied by other companies.
3) Add the separation of resin from metal frames
- Use spherical retention beads with a particle size of 100 to 200 µm on the facing surface of the metal frame wax mold.
- Sandblast the facing surface of the metal frame with alumina particles, 50 µm diameter.
- If the metal frame is a precious metal alloy, apply a metal adhesive primer (eg. ALLOY PRIMER) or use a dental tin plating unit to coat the metal frame.
- To securely harden the Opaque resin and assure strong adhesion of the Opaque resin to the metal, first apply Opaque Primer to the adhering surface of the metal. Apply the Opaque resin after being sure the Opaque Primer has dried completely.
- The standard procedure is to apply the Opaque twice. Be sure the first layer of the Opaque has hardened completely before applying Opaque Primer on the first layer of the Opaque. Dry the second Primer coat, and then apply a second layer of the Opaque.
- Use Opaque Modifier for color adjustment of the Body Opaque.
4) Avoid separation of the Body resin from the EG Fiber
- Sandblast the EG Fiber with alumina particles, 50 µm diameter.
- Do not touch the EG Fiber with your bare hands or fingers when sandblasting or applying the primer, the fiber surface will be contaminated, preventing the EG Fiber from adhering to the Body.
5) Avoid roughness in the oral cavity
- When heat curing the pastes coats, be sure to use the heating temperature and length of time specified in these Instructions for Use.
- To minimize the formation of unpolymerized resin, which leads to a rough surface, heat curing should be performed in the 100°C-110°C/212°F-230°F range for 15 minutes.
- Before final light curing, apply the Air Barrier Paste to the restoration surface and keep it on until the heat curing is over.
- Shaping, finishing and polishing should be performed after heat curing. If heat curing is performed after shaping, finishing or polishing, insufficient surface hardness could result, thus leading to premature roughening of the restoration while in the oral cavity.
- Poor initial smoothness will cause premature surface roughness after the restoration is cemented in the oral cavity. Observe the following conditions:
1. Remove rough scores and scratches with a silicone point before polishing.
2. When polishing the restoration, use the polishing agent, brush and felt wheel supplied with the product.
3. The restoration should be, and can be, polished until the surface is as smooth as glazed porcelain.




7) Avoid poor polymerization
- The Opacous Dentin should be built up in thin layers and each layer light cured for the specified length of time; otherwise poor polymerization will result.
- After curing the EG Fiber using EG Core instrument, remove the EG Core instrument and then light cure the EG Fiber again to ensure optimum polymerization.
- In very deep cavity preparations, the inner surface at the floor of the cavity may not cure sufficiently. If the cavity is deeper than the polymerization depth indicated, build up the paste layers incrementally, at light cure each installation.

8) Avoid poor separation from plaster models
- Use the separating agent supplied with the product kit; otherwise the restoration may not release properly or the plaster model may break.
- Be careful not to let the cured resin extend over the margins; otherwise the restoration may not release properly or the plaster model may break during release. Scrape the cured resin from the margins before releasing the restoration from the model.

9) Avoid polymerization of pastes in ambient light
- Do not mix different types of Body resin. If they are used under strong light, such as by a window or under a laboratory light. Use the materials in more moderate light conditions.
10) Avoid entrapment air bubbles
- Apply Modeling Liquid on layered surfaces whenever adding layering Body resin.
- Do not mix different types of Body resin.
- Be careful to prevent air bubbles from being trapped between the EG Fiber and the Body resin when placing the Body resins on the EG Fiber.

11) Other points to avoid
- Do not use the same brush to apply Opaque Primer and Opaque resin; the Opaque resin will cure on the brush before application when the same brush is used. Brushes used to apply Opaque resin and Opaque Primer should be washed in ethyl alcohol or a self-curing resin monomer liquid after use.
- Be careful not to burn yourself when softening EG Core in hot water.
12) Storage precautions
- Do not leave the caps off the Opaque Primer, Add-On Primer and CR Sep. After use, the cap should be replaced as soon as possible to prevent premature setting or contamination.
- The EG Fiber should be put into an aluminum bag very soon after use to prevent it from hardening.
- The Opaque Primer, Add-On Primer and CR Sep are flammable. Do not use or store them near an open flame.
- The product should not be stored in direct sunlight or near an open flame. It should be stored at room temperature (2°C-25°C/39°F-77°F).
13) Expiration date
The product should be used by the expiration date* indicated on the package. [The expiration date is established by Kuraray Medical's own certification (from Kuraray Medical research data). * (Example: 2005-11 means that the expiration date is Nov 2005.)

V. COMPONENTS AND CHEMICAL INGREDIENTS
V-1. Components
- System components of ESTENIA™ C&B: They are shown each on package.
- Single items: The following items are available as a single component of ESTENIA™ C&B.
Each shade of Body resin and Opaque resin, EG Fiber (For anterior teeth, For posterior teeth), EG Flow, Opaque Primer, Modeling Liquid, Add-On Primer, Jacket Separate Kit (Jacket Spacer and Margin Sep), Jacket Spacer, Margin Sep, CR Sep, Air Barrier Paste, Polishing set (Polishing compound, Polishing brush, Felt wheel, EG Core, Instruments (No.1, No.2, No.3 or set), Mixing plate (#954), Small brush (#901), Small brush (#932), Small brush (#933)
V-2. Chemical ingredients
- Body Resin: Monomer (Polyurethane methacrylonomer and methacrylic acid series monomer), Filler (Surface treated glass powder and surface treated aluminum micro filler), Photocuring catalyst, Colorant and others
- Opaque Resin: Monomer (Bis-GMA and other methacrylic acid series monomer), Filler (Surface treated quartz powder, Prepolymerized resin filler and others), Photocuring catalyst, Colorant and others
- Opaque Primer: Methacrylic acid series monomer, MDP, Solvent and others
- Modeling Liquid: Methacrylic acid series monomer, Polyurethane methacrylonomer, MDP, Photocuring catalyst and others
- Add-On Primer: Silane coupling agent and ethanol
- EG Fiber: Monomer (Polyurethane methacrylonomer, TEGDMA and other methacrylic acid series monomer), Surface treated glass filler, Surface treated micro filler, Photocuring catalyst, Colorant and others
V-3. Basic laboratory procedures
1) Fabrication of jacket crowns
1) Fabrication of a working model
Fabricate a working model in the usual manner, and then trim the model.
2) Application of the separating agent
1) Applying and drying Jacket Spacer: To make it easy to release the crown from the mold, apply Jacket Spacer uniformly to all the area, except the margins, with a brush, and then dry the applied Jacket Spacer by leaving it or by blowing air gently over the coating.
2) Applying and drying Margin Sep: Apply Margin Sep on the margins of the abutment and dry by leaving it or by blowing air gently over it.
3) Forming and preliminary polymerization of Opacous Dentin (when abutment is composite resin)
When the abutment is composite resin, use Opacous Dentin for the substrate color.
4) Application and light curing of Opaque resin (when abutment is metal)
When the abutment is metal, use Opaque resin.
1) Applying and drying Opaque Primer: Apply Opaque Primer on the abutment tooth and evaporate the volatile content of the primer by mildly blowing with air or leaving it for about 30 seconds.
2) Applying and light curing Opaque resin: Apply Body Opaque paste of the shade specified on the abutment surface thinly and light cure the resin for the specified length of time using a laboratory-type light curing unit. When you want to enhance the color in the cervical area, apply Cervical Opaque in the cervical area of the abutment. Opaque Modifier may be used to adjust the Body Opaque color by mixing with Body Opaque paste to apply it alone.
3) Applying and light curing the 2nd layer of Opaque Primer and Opaque resin: Apply Opaque Primer on the cured 1st layer of Body Opaque and dry it. Then, apply another layer of Body Opaque and light cure it in the same manner as for the 1st layer. Repeat this process until the spacer color is masked out.
5) Forming and preliminary polymerization of Body resin
1) Forming and preliminarily light curing Cervical: As necessary, form the paste to an area from the margin to 1/4 - 1/5 of the crown length to the cervical area, making the paste thickness gradually thinner. Then, light cure the resin preliminarily for the specified length of time.
2) Applying Modeling Liquid: Apply a thin coat of Modeling Liquid on the surface of the cured resin for better wetting and improved adaptation when layering another layer of paste. Do not apply Modeling Liquid too much; otherwise Dentin paste might stick to the cured resin. Wipe away any excess liquid with tissues.
3) Forming and preliminarily light curing Dentin: To produce the basic color of the crown, form a generous amount of Dentin paste and light cure it preliminarily for the specified length of time.
4) Forming and preliminarily light curing Transparent: Apply a thin coat of Transparent paste to 1/4 to 1/5 of the crown length to the occlusal surface, and then light cure the paste preliminarily for the specified length of time.
5) Forming and preliminarily light curing Enamel: Form Enamel paste to 1/3 of the crown length to the occlusal surface and then light cure the paste preliminarily for the specified length of time.
6) Application of Air Barrier paste
After forming Enamel paste in the contact area, apply Air Barrier paste on the surface layer of the crown to reduce the generation of unpolymerized resin.
7) Final light curing
Light cure the pastes for the specified length of time.
8) Heat curing
After releasing the crown from the mold, heat cure the restoration for 15 minutes at a temperature of 100°C-110°C/212°F-230°F.
9) Shaping and finishing
Shape using a carborundum point or other instruments, then remove scratches and grooves thoroughly from the crown surface using a silicone point.
10) Glazing
Glaze the restoration using the polishing set supplied with the kit.
1) Polishing with brush: Apply the polishing compound to the brush and polish the occlusal surface and areas around it with the brush.
2) Finishing with brush: Apply the polishing compound to the felt wheel and glaze the entire crown using the felt wheel.
11) Completion
Sandblast the inner surface of the crown using alumina particles, 30-50 µm diameter, at a pressure of 1 to 2 kg/cm².
2) Fabrication of inlays and onlays
1) Fabrication of a working model
If there is an undercut in the cavity, block it out using plaster or wax.
2) Application and drying of the resin separating agent
Apply CR Sep on trays around the cavity, adjacent teeth or opposing tooth.
3) Forming and preliminary polymerization of Body resin
1) Application and preliminary polymerization of Dentin: If the cavity preparation has a proximal area, form Dentin paste in that proximal area and light cure the paste preliminarily for the specified length of time.
2) Application and preliminary polymerization of Cervical Transparent: Form Cervical Transparent paste on the occlusal surface of the cavity and then light cure the paste preliminarily for the specified length of time.
3) Forming and preliminary polymerization of Enamel: Form Enamel paste on the occlusal surface and light cure the paste preliminarily for the specified length of time.

4) Application of Air Barrier paste and final light curing
Apply Air Barrier paste on the surface layer of the inlay, then light cure the pastes for the specified length of time.
5) Heat curing
After releasing the inlay from the mold, heat cure the restoration for 15 minutes at a temperature of 100°C-110°C/212°F-230°F.
6) Shaping, finishing, glazing and completion
As in making a jacket crown, shape, finish and glaze the restoration. Then remove the separating agent and sandblast the inner surface of the inlay for pretreatment for adhesion.
3) Fabrication of bridges with frameworks
Be careful not to expose EG Fiber on the restoration's surface when fabricating a framework.
1) Fabrication of a working model
To provide a uniform space between the framework and the mucosal surface of the pontic, keep the pontic about 1 mm above the mucosal level.
2) Making of the frame fabrication core
1) Waxing up: Using the sprue wax (about 2 mm diameter for the anterior region, about 2.5 - 3 mm diameter for the posterior region), reproduce the framework contour. The framework contour should be designed so that it is located under the pontic as much as possible. Check the occlusal relationship to be sure sufficient gap is attained.
2) Blocking out with silicone putty: To prevent the framework from displacing due to placement of EG Core, block out the area under the wax frame and the axial surface of the abutment with silicone putty.
3) Softening EG Core: Cut EG Core into a proper length and immerse it in a hot water bath (about 80°C/176°F) for about 3 minutes to soften the Core material. Then place it onto the wax frame of the abutment.
4) Completion of the frame fabrication core: After curing, remove EG Core from the model.
3) Fabrication of the framework
1) Applying the separating agent: Apply CR Sep on the plaster model.
2) Preparing EG Fiber: Cut EG Fiber into a size long enough to cover the entire occlusal surface of the bridge.
3) Placing EG Fiber into the inner surface of the core: Place EG Fiber into the inner surface of EG Core or on the model.
4) Placing the core under pressure: Return EG Core onto the model. After that, light cure the core preliminarily for 60 seconds using the laboratory-type light curing unit.
5) Light curing the framework: Remove the core from the model and light cure it for the specified length of time to cure EG Fiber.
6) Adjusting the framework: Remove EG Fiber from the model and trim using a carborundum wheel. Return EG Fiber onto the model and check the occlusal surface for proper clearance.
4) Application of the resin separating agent
Apply Jacket Spacer uniformly to all the area of the abutment except the margins with a brush and Margin Sep on the margins, and dry it by leaving it or by blowing air gently.
5) Drying of the framework
1) Sandblasting: Sandblast the entire framework with alumina particles.
2) Bonding treatment: Apply Add-On Primer on the framework surface, After drying it, apply Modeling Liquid.
3) Forming Opacous Dentin: Apply a thin coat of Opacous Dentin onto the basal surface of the pontic and light cure the paste for the specified length of time. Form the paste thinly on the abutment. If the abutment is metal, it is advisable to apply Opaque Primer and Opaque resin, and light cure the pastes for the specified length of time in advance.
4) Placing the framework: Form Dentin or Opacous Dentin on the pontic area and place the framework on the abutment area on the model.
5) Applying EG Flow: Apply EG Flow on the connection areas of the frame and the joint part of the framework and paste.
6) Light curing: Light cure the paste for the specified length of time to secure the framework to the abutment.
6) Forming Body resin
As in making a jacket crown, form Dentin, Transparent and Enamel.
7) Applying Air Barrier paste and final light curing
Apply Air Barrier paste on the surface layer of the bridge, then light cure the pastes for the specified length of time.
8) Heat curing
After releasing the bridge from the mold, heat cure the restoration for 15 minutes at a temperature of 100°C-110°C/212°F-230°F.
9) Shaping, finishing, glazing and completion
As in making a jacket crown, shape, finish and glaze the restoration. Then sandblast the inner surface of the crown for pretreatment for adhesion.
4) Fabrication of facing cast crown and facing cast bridge
1) Fabrication of a working model
Fabricate a working model in the usual manner.
2) Fabrication of metal frame
Fabricate a metal frame by waxing up and casting in the usual manner.
Apply retention beads with a diameter of 100 - 200 µm to the facing surface.
3) Pretreatment of metal frame
1) Sandblasting: Sandblast the metal frame with alumina particles, 50 µm diameter. After that, clean the surface with ultrasound for 2 minutes and dry.
2) Pretreatment for metal adhesive: If a precious metal alloy is used, apply the metal adhesive primer (eg. ALLOY PRIMER) to the facing surface of the metal frame before drying.
3) Apply and light curing Opaque resin
1) Apply and dry Opaque Primer: Apply Opaque Primer to the facing surface and dry the primer by mildly blowing with air or leaving it for about 30 seconds to evaporate the volatile content of the primer.
2) Apply and light cure Opaque resin: Apply Body Opaque paste of the shade specified on the abutment surface thinly and light cure the resin for the specified length of time using a laboratory-type light curing unit.
3) Apply and light cure the 2nd layer of Opaque Primer and Opaque resin: Apply Opaque Primer on the cured 1st layer of Body Opaque and dry it. Then, apply another layer of Body Opaque and light cure it in the same manner as for the 1st layer. Repeat this process until the spacer color is masked out.
4) Completion of the crown
As in making a jacket crown, after forming of Body resin, light cure, heat cure, shape, finish, and glaze to complete the crown.

WARANTY
KURARAY MEDICAL INC. will replace any product that is proved to be defective. KURARAY MEDICAL INC. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.
(NOTE)
ESTENIA is a trademark of KURARAY CO., LTD.

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1) Colori

Color chart table with columns A1-A4, B1-B4, C1-C4 and rows for Body Opaque, Cervical Opaque, Opaque Modifier, Opaque Dentin, Dentin, Cervical, Enamel, Transparent, Cervical Transparent, Cervical Dentin Effect, Enamel Effect, Transparent Effect, Gingival.

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2) Colori

Color chart table with columns C2-C4, D2-D4, NWO, NWO.5, NP1.5, NP2.5 and rows for Body Opaque, Cervical Opaque, Opaque Modifier, Opaque Dentin, Dentin, Cervical, Enamel, Transparent, Cervical Transparent, Cervical Dentin Effect, Enamel Effect, Transparent Effect, Gingival.

3) Colori

Color chart table with columns C2-C4, D2-D4, NWO, NWO.5, NP1.5, NP2.5 and rows for Body Opaque, Cervical Opaque, Opaque Modifier, Opaque Dentin, Dentin, Cervical, Enamel, Transparent, Cervical Transparent, Cervical Dentin Effect, Enamel Effect, Transparent Effect, Gingival.

4) Colori

Color chart table with columns C2-C4, D2-D4, NWO, NWO.5, NP1.5, NP2.5 and rows for Body Opaque, Cervical Opaque, Opaque Modifier, Opaque Dentin, Dentin, Cervical, Enamel, Transparent, Cervical Transparent, Cervical Dentin Effect, Enamel Effect, Transparent Effect, Gingival.

5) Colori

Color chart table with columns C2-C4, D2-D4, NWO, NWO.5, NP1.5, NP2.5 and rows for Body Opaque, Cervical Opaque, Opaque Modifier, Opaque Dentin, Dentin, Cervical, Enamel, Transparent, Cervical Transparent, Cervical Dentin Effect, Enamel Effect, Transparent Effect, Gingival.

6) Colori

Color chart table with columns C2-C4, D2-D4, NWO, NWO.5, NP1.5, NP2.5 and rows for Body Opaque, Cervical Opaque, Opaque Modifier, Opaque Dentin, Dentin, Cervical, Enamel, Transparent, Cervical Transparent, Cervical Dentin Effect, Enamel Effect, Transparent Effect, Gingival.

7) Colori

Color chart table with columns C2-C4, D2-D4, NWO, NWO.5, NP1.5, NP2.5 and rows for Body Opaque, Cervical Opaque, Opaque Modifier, Opaque Dentin, Dentin, Cervical, Enamel, Transparent, Cervical Transparent, Cervical Dentin Effect, Enamel Effect, Transparent Effect, Gingival.

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