

CZR PRESS PRESS CERAMIC for Zirconia Framework Restorations

CZR PRESS SHADE BASE STAIN KIT

CE 0120

CZR PRESS SHADE BASE STAIN KIT is a kind of the stain to make the basic shade color of zirconia framework before pressing CZR PRESS INGOT. The basic shade color can be recreated thinly, transparently and naturally with the simple stain technique.

BASIC SHADE	3g each	20 colors
IS LIQUID	10ml	1 bottle
SS COLOR GUIDE		1 piece
TECHNICAL INSTRUCTIONS		1 copy



CZR PRESS SS MODIFIER KIT

CE 0120

CZR PRESS SS MODIFIER KIT is assortment of all available SS Modifiers. It can be mixed with CZR PRESS Shade Base Stain to customize the shade or may be applied alone for minor modifications.

CV ORANGE	A+	3g each
INCISAL BLUE 1	B+	
INCISAL BLUE 2	C+	
GRAY	D+	
SALMON PINK	FLUORO	
EARTH BROWN		
IS LIQUID	10ml	1 bottle
SS MODIFIER COLOR GUIDE		1 piece
TECHNICAL INSTRUCTIONS		1 copy



CZR PRESS INGOT-L KIT / INGOT-H KIT

CE 0120

CZR PRESS INGOT-L KIT is a low translucent type of press ceramic for zirconia framework. It is developed to recreate the margin area precisely by pressing INGOT-L. INGOT-L is a press ceramic using with the layering technique. After INGOT-L is pressed, the crown can be made by using CZR TRANSLUCENT PORCELAIN or CZR LUSTER PORCELAIN. Can make a single anterior crown without frameworks using with CZR PRESS LF porcelain.

BASIC SHADE	10g each	20 colors
PLUNGER		1 set
TECHNICAL INSTRUCTIONS		1 copy

5g Ingot is available as individual item.



CZR PRESS INGOT-L KIT

CZR PRESS GLAZE

CE 0120

This CZR PRESS GLAZE POWDER is necessary to produce glossiness of the crown in the H-ingots staining technique. Be sure to use this glaze powder and refer to our CZR PRESS technical instruction for the glazing temperature. This glaze powder is available individually. **The use of IS Liquid is highly recommended.**

CZR PRESS GLAZE	10g
-----------------	-----



CZR PRESS PRESS CERAMIC for Zirconia Framework Restorations

CZR PRESS TRIAL-L KIT

CE 0120

CZR PRESS TRIAL-L KIT is a necessary kit for using INGOT-L, including SHADE BASE STAIN, INGOT-L, CZR PORCELAIN, and PLUNGER. It contains 5 shade colors, which are A1, A2, A3, NW0, and NW0.5. The beginning user of CZR PRESS should also prepare for press investments, RING/FORMER, IS LIQUID, and CZR FORMING LIQUID.

L-Ingot

for use when utilizing the "Layering Method" & "LF Layering Method"

INGOT-L	10g each	5 colors
SB STAIN	3g each	5 colors
CZR PORCELAIN POWDER	10g each	8 colors
PLUNGER		1 set
SB STAIN COLOR GUIDE		1 piece
CZR PORCELAIN INSTRUCTION		1 copy
TECHNICAL INSTRUCTIONS		1 copy



CZR PRESS TRIAL-H KIT

CE 0120

CZR PRESS TRIAL-H KIT is a necessary kit for using INGOT-H, including SHADE BASE STAIN, INGOT-H, CZR PRESS GLAZE, CZR ES (A+), and PLUNGER. It contains 5 shade colors, which are A1, A2, A3, NW0, and NW0.5. The beginning user of CZR PRESS should also prepare for press investments, RING/FORMER, IS LIQUID, and ES LIQUID.

H-Ingot

for use when utilizing the "Staining Method" & "LF Layering Method"

INGOT-H	10g each	5 colors
SB STAIN	3g each	5 colors
CZR ES A+	3g	1 color
CZR PRESS GLAZE	10g	1 color
PLUNGER		1 set
SB STAIN COLOR GUIDE		1 piece
CZR PORCELAIN INSTRUCTION		1 copy
TECHNICAL INSTRUCTIONS		1 copy



CZR PRESS ESTHETIC WHITE KIT

CE 0120

CZR PRESS ESTHETIC WHITE KIT is developed to create whiter shades rather than NW shades with Press technique. It is not only pure bleaching shades but also natural and beautiful white shades.

SHADE BASE STAIN	3g each	SS White, SSB1
ESTHETIC WHITE INGOT	10g each	EW00, EW0, EW, EWY
EXTERNAL STAIN	3g	B+
GLAZE	10g	CZR PRESS Glaze
TECHNICAL INSTRUCTION		1 copy



CZR PRESS LF BASIC KIT

CE 0120

CZR PRESS LF BASIC KIT is basic assortment for layering for stand-alone CZR PRESS crowns. It can be used with normal CZR layering also. 7 LF BODY shades are very useful for minor correction of pressed body shape.

BODY	10g each	HA1B, HA3B, HA4B HB2B, HC2B, HD2B HNW0B
ENAMEL		E1, E2, E3
TRANSLUCENT		Tx, To, T1, T2
MODIFIER		Mamelon1~2
ADD-ON		AD-T, AD-B
MARGIN		MRP



CZR PRESS PRESS CERAMIC for Zirconia Framework Restorations

CZR PRESS LF LUSTER CCV KIT

CE 0120

CZR PRESS LF LUSTER/CCV KIT is assortment of all available CZR PRESS LF LUSTER and CCV porcelain. LUSTER can produce beautiful surface texture with natural opalescent. CCV is very effective when applying on cervical area.

LTo	LT ₁	T Blue	AQUA BLUE 1	10g each
AQUA BLUE 2	LT NATURAL	LT SUPER GRAY	CREAMY ENAMEL	
SUN BRIGHT	INCISAL AUREOLA	CREAMY WHITE	LY YELLOW	
CCV-1	CCV-2	CCV-3	CCV-4	



CZR PRESS LF ESTHETIC WHITE LAYERING KIT

CE 0120

CZR PRESS LF ESTHETIC WHITE LAYERING KIT is developed to create whiter shades rather than NW shades. It is not only pure bleaching shades but also natural and beautiful white shades.

SHADE BASE STAIN	3g each	SB White, SBB ₁
BODY	10g each	EW ₀₀ B, EW ₀ B, EW _B , EWYB
EXTERNAL STAIN	3g	B ⁺
GLAZE	10g	CZR PRESS Glaze
TECHNICAL INSTRUCTION	1 copy	



CZR PRESS LF TISSUE KIT

CE 0120

CZR PRESS LF TISSUE KIT is complete kit of 7 Noritake Tissue colors. All colors are very useful for implant work.

TISSUE	10g each	Tissue1~7
--------	----------	-----------



CZR PRESS LF INTERNAL STAIN COMPLETE KIT

CE 0120

CZR PRESS LF INTERNAL STAIN KIT is a set of Low fusing stains for CZR PRESS LF porcelain. It is necessary to use this stain before CZR PRESS LF layering. Baking temperature of this stain is 800°C.

WHITE	MAMELON ORANGE 1	CERVICAL 2	REDDISH BROWN	GRAY	B ⁺	3g each
INCISAL BLUE 1	MAMELON ORANGE 2	CERVICAL 3	SALMON PINK	FLUORO	C ⁺	
INCISAL BLUE 2	CERVICAL 1	EARTH BROWN	RED	A ⁺	D ⁺	
BRIGHT (DILUTION)					3g	1 color
IS LIQUID					10ml	1 bottle
IS COLOR GUIDE						1 piece
TECHNICAL INSTRUCTIONS						1 copy



CZR PRESS LF EXTERNAL STAIN KIT

CE 0120

CZR PRESS LF EXTERNAL STAIN KIT is a set of Low fusing stains for CZR PRESS LF porcelain. It is necessary to use this stain after CZR PRESS LF layering. Baking temperature of this stain is 840°C.

PURE WHITE	BLUE	YELLOW	CERVICAL 1	EARTH BROWN	PINK	B ⁺	3g each
GRAY	GREEN 1	ORANGE 1	CERVICAL 2	REDDISH BROWN	RED	C ⁺	
BLACK	GREEN 2	ORANGE 2	CERVICAL 3	SALMON PINK	A ⁺	D ⁺	
GLAZE					10g	1 color	
ES LIQUID					10ml	1 bottle	
ES COLOR GUIDE						1 piece	
TECHNICAL INSTRUCTIONS						1 copy	



RING / FORMER SET

NORITAKE RING / FORMER SET is used for investing waxed-up crowns. 3 sizes are available, 100g, 200g and 300g.

LF LIQUID

CZR PRESS LF porcelain should EXCLUSIVELY be mixed with LF LIQUID.

NORITAKE DISPO PLUNGER

Noritake DISPO PLUNGER is disposable plunger set. 2 kinds of size are available for each 2g and 5g CZR PRESS INGOT.