

VITA 3D Master* Renk Skalası ile Kuraray Noritake EX-3 Dönüşüm Tablosu

RENK	OPAK	DENTİN	MAKYAJ
1M1	A ₁	A ₁	
1M2	A ₁	A ₁ : A ₂ = 1:1	
2L1.5	A ₂	A ₂ : A ₃ = 1:1	
2L2.5	A ₃	A ₃	Sarı
2M1	C ₁	A ₁	
2M2	A ₂	A ₃	
2M3	B ₃	B ₃	
2R1.5	D ₂	A ₂	
2R2.5	A ₃	A ₃	
3L1.5	D ₃	D ₃	Kahverengi
3L2.5	A ₃	A ₃	Kahverengi
3M1	C ₁	C ₁	Kahverengi
3M2	A ₃	A ₃	Kahverengi
3M3	B ₄	B ₄	
3R1.5	D ₃	D ₂	
3R2.5	A _{3.5}	A ₃ : A _{3.5} = 1:1	
4L1.5	C ₂	C ₂ : Mod-Brown = 10:1	
4L2.5	A ₄	A ₄	Turuncu
4M1	C ₂	C ₁ : Mod-Brown = 8:1	
4M2	A _{3.5}	A ₃ : A _{3.5} = 1:3	
4M3	A ₄	C ₂ :Mod-Orange:Mod-Brown = 100 : 8 : 1	
4R1.5	A ₄	C ₃ : Mod-Brown = 8:1	
4R2.5	A ₄	A ₄	Kahverengi
5M1	C ₃	C ₃ : Mod-Brown = 8:1	
5M2	A ₄	C ₄ : Mod-Brown = 10:1	
5M3	B ₄	A ₄ : Mod-Orange = 8:1	