

Super Porcelain EX-3 PRESS PRESS CERAMIC for Metal Framework Restorations

NORITAKE SUPER PORCELAIN

EX-3 PRESS PASTE OPAQUE BASIC KIT

CE 0120

EX-3 PRESS PASTE OPAQUE BASIC KIT is specially designed for EX-3 PRESS INGOT. It can make the basic shade color of metal framework with a thinner coat before pressing ingot.

PPnA1	PPnB1	PPnC2	PPnD4	6g each
PPnA2	PPnB2	PPnC3	PPNW0	
PPnA3	PPnB3	PPnC4	PPNW0.5	
PPnA3.5	PPnB4	PPnD2	PPNP1.5	
PPnA4	PPnC1	PPnD3	PPNP2.5	
PPEW0	PPEW			
TECHNICAL INSTRUCTIONS				1 copy
PASTE OPAQUE MODIFIER BRUSH				1 piece



NORITAKE SUPER PORCELAIN

EX-3 PRESS PASTE OPAQUE MODIFIER KIT

CE 0120

NORITAKE EX-3 PRESS PASTE OPAQUE MODIFIER KIT contains 8 color modifiers which may be incrementally incorporated into EX-3 PRESS Basic Paste Opaques to create custom shades. They may also be applied individually for minor color adjustments & subtle shade variances.

PP WHITE	PP REDDISH BROWN	3g each
PP GRAY	PP PINK	
PP ORANGE	PP BLUE	
PP EARTH BROWN	PP YELLOW	
TECHNICAL INSTRUCTIONS		1 copy
PASTE OPAQUE MODIFIER BRUSH		1 piece



NORITAKE SUPER PORCELAIN

EX-3 PRESS INGOT KIT L / H

CE 0120

NORITAKE EX-3 PRESS INGOT KIT L contains low translucency pressable ceramic ingots. It was developed to recreate the margin area precisely by pressing INGOT L. The L ingot is indicated for use in the layering technique. After ingot L is pressed, the crown can be made with EX-3 PRESS LF Enamel, LF Translucent or Luster porcelain, LF Internal stain and LF External stain.

INGOT KIT L

HnA1	HnB1	HnC2	HnD4
HnA2	HnB2	HnC3	HNW0
HnA3	HnB3	HnC4	HNW0.5
HnA3.5	HnB4	HnD2	HNP1.5
HnA4	HnC1	HnD3	HNP2.5
TECHNICAL INSTRUCTIONS			1 copy
EX-3 PRESS DISPO PLUNGER			3 pieces

INGOT KIT H

LnA1	LnA2	LnA3	LnA3.5
LnA4	LnB1	LnB2	LnB3
LnB4	LnC1	LnC2	LnC3
LnC4	LnD2	LnD3	LnD4
LNW0	LNW0.5	LNP1.5	LNP2.5
TECHNICAL INSTRUCTIONS			1 copy
EX-3 PRESS DISPO PLUNGER			3 pieces



EX-3 PRESS INGOT KIT L

NORITAKE SUPER PORCELAIN

EX-3 PRESS ESTHETIC WHITE KIT

CE 0120

EX-3 PRESS Esthetic White Kit is for creating whiter shades rather than NW shades, which is not only pure bleaching shades but also natural and beautiful white shades.

PPEW0	PPEW	EW00	EW0
EW	EWY	LF Silky E1	LF Silky E2
LF ELT1	LF ELT2	LF ELT3	
TECHNICAL INSTRUCTIONS			1 copy
EX-3 PRESS DISPO PLUNGER			3 pieces



Super Porcelain EX-3 PRESS PRESS CERAMIC for Metal Framework Restorations

NORITAKE SUPER PORCELAIN

EX-3 PRESS TRIAL KIT L1 / L2

CE 0120

EX-3 PRESS TRIAL KIT L1 & L2 are an introductory kit which includes EX-3 PRESS Opaque, L ingots, EX-3 PRESS LF Enamel, LF Translucent, LF Luster, AD-T, AD-B and MRP Porcelains. TRIAL KIT L1 contains the following 5 shades: NW₀, NW_{0.5}, nA₁, nA₂ & nA₃. And TRIAL KIT L2 contains the following 5 shades: nA₂, nA₃, nA_{3.5}, nB₂ & nC₂. For EX-3 PRESS procedures, Noritake PRESS Investment, Ring / former, Paste Opaque Liquid & ES Liquid are required.

TRIAL KIT L1

PPnA ₁	PPnA ₂	PPnA ₃	PPNW ₀
PPNW _{0.5}	LnA ₁	LnA ₂	LnA ₃
LNW ₀	LNW _{0.5}	LF Silky E ₂	LF E ₂
LF E ₃	LF LT ₁	LF CE	LF Tx
LF AD-T	LF MRP		
TECHNICAL INSTRUCTIONS			1 copy
EX-3 PRESS DISPO PLUNGER			1 piece

TRIAL KIT L2

PPnA ₂	PPnA ₃	PPnA _{3.5}	PPnB ₂
PPnC ₂	LnA ₂	LnA ₃	LnA _{3.5}
LnB ₂	LnC ₂	LF E ₂	LF E ₃
LF LT ₁	LF CE	LF Tx	LF AD-T
LF AD-B	LF MRP		
TECHNICAL INSTRUCTIONS			1 copy
EX-3 PRESS DISPO PLUNGER			3 piece



EX-3 PRESS TRIAL KIT L1

NORITAKE SUPER PORCELAIN

EX-3 PRESS TRIAL KIT H1 / H2

CE 0120

EX-3 PRESS TRIAL KIT H1 & H2 are an introductory kit which includes EX-3 PRESS Opaque, H ingots, EX-3 PRESS LF External Stain A⁺ & EX-3 PRESS LF Glaze Powder. TRIAL KIT H1 contains the following 5 shades: NW₀, NW_{0.5}, nA₁, nA₂ & nA₃. And TRIAL KIT H2 contains the following 5 shades : nA₂, nA₃, nA_{3.5}, nB₂ & nC₂. The beginner EX-3 PRESS user should also use Noritake PRESS Investment, Ring / former, Paste Opaque Liquid & ES Liquid.

TRIAL KIT H1

HnA ₁	HNW ₀	PPnA ₂	PPNW _{0.5}
HnA ₂	HNW _{0.5}	PPnA ₃	A ⁺
HnA ₃	PPA ₁	PPNW ₀	GLAZE
TECHNICAL INSTRUCTIONS			1 copy
EX-3 PRESS DISPO PLUNGER			3 pieces

TRIAL KIT H2

HnA ₂	HnC ₂	PPnB ₂	C ⁺
HnA ₃	PPnA ₂	PPnC ₂	GLAZE
HnA _{3.5}	PPn ₃	A ⁺	
HnB ₂	PPn _{3.5}	B ⁺	
TECHNICAL INSTRUCTIONS			1 copy
EX-3 PRESS DISPO PLUNGER			3 pieces



EX-3 PRESS TRIAL KIT H1

NORITAKE SUPER PORCELAIN

EX-3 PRESS LF BODY KIT

CE 0120

EX-3 PRESS LF body can be used for repairing and for additional touch-up on pressed area. This Body Kit is composed of All n Colors, EW, NW and NP shades.

LF nA ₁ B	LF nA ₂ B	LF nA ₃ B	LF nA _{3.5} B	10g each
LF nA ₄ B	LF nB ₁ B	LF nB ₂ B	LF nB ₃ B	
LF nB ₄ B	LF nC ₁ B	LF nC ₂ B	LF nC ₃ B	
LF nC ₄ B	LF nD ₂ B	LF nD ₃ B	LF nD ₄ B	
LF EW ₀ B	LF EW ₀ B	LF EW ₁ B	LF EW ₂ B	
LF NW ₀ B	LF NW _{0.5} B	LF NP _{1.5} B	LF NP _{2.5} B	



NORITAKE SUPER PORCELAIN

EX-3 PRESS LF COMPLETE KIT

CE 0120

EX-3 PRESS LF is low fusing porcelain to be built up not only to make an enamel layer but also cervical layer after pressing EX-3 press ingot. This kit contains LF Enamel, Translucent, Luster, Clear Cervical, Modifier, Add-on, MRP porcelains. By using these porcelains with pressable ceramic, you can make more natural crown.

LF E ₁	LF E ₂	LF E ₃	LF Tx	LF T ₀	10g each
LF T ₁	LF T ₂	LF LT ₀	LF LT ₁	LF T-BLUE	
LF AQUA BLUE 1	LF AQUA BLUE 2	LF LT NATURAL	LF LT SUPER GRAY	LF CREAMY ENAMEL	
LF SUN BRIGHT	LF INCISAL AUREOLA	LF CREAMY WHIT	LF LT YELLOW	LF CCV-1	
LF CCV-2	LF CCV-3	LF CCV-4	LF Mamelon 1	LF Mamelon 2	
LF AD-T	LF MRP				
TECHNICAL INSTRUCTIONS					



Super Porcelain EX-3 PRESS PRESS CERAMIC for Metal Framework Restorations

NORITAKE SUPER PORCELAIN

EX-3 PRESS LF ESTHETIC WHITE LAYERING KIT

CE 0120

EX-3 PRESS LF ESTHETIC WHITE LAYERING KIT is assortment of Esthetic White shades of EX-3 PRESS LF porcelains. It can crate whiter shades rather than NW shades.

PPEW ₀	PPEW	LF EW ₀₀ B	LF EW ₀ B
LF EW _B	LF EW _{YB}	LF Silky E ₁	LF Silky E ₂
LF ELT ₁	LF ELT ₂	LF ELT ₃	
PO BRUSH			1 piece
TECHNICAL INSTRUCTIONS			1 copy



NORITAKE SUPER PORCELAIN

EX-3 PRESS LF TISSUE KIT

CE 0120

EX-3 PRESS LF TISSUE KIT is complete kit of 7 Noritake Tissue colors. All colors are very useful for implant work.

LF TISSUE 1	LF TISSUE 2	LF TISSUE 3	LF TISSUE 4
LF TISSUE 5	LF TISSUE 6	LF TISSUE 7	
TECHNICAL INSTRUCTIONS			1 copy



EX-3 PRESS LF INTERNAL STAIN COMPLETE KIT

CE 0120

EX-3 PRESS LF INTERNAL STAIN KIT is Low fusing stain for EX-3 PRESS LF porcelain. It is necessary to use this stain before EX-3 PRESS LF layering

LF WHITE	LF CERVICAL 1	LF SALMON PINK	LF B+	3g each
LF INCISAL BLUE 1	LF CERVICAL 2	LF RED	LF C+	
LF INCISAL BLUE 2	LF CERVICAL 3	LF GRAY	LF D+	
LF MAMELON ORANGE 1	LF EARTH BROWN	LF FLUORO	LF BRIGHT	
LF MAMELON ORANGE 2	LF REDDISH BROWN	LF A+		
IS LIQUID			10ml	1 bottle
IS COLOR GUIDE				1 piece
TECHNICAL INSTRUCTIONS				1 copy



EX-3 PRESS LF EXTERNAL STAIN KIT

CE 0120

EX-3 PRESS LF EXTERNAL STAIN KIT is a set of Low fusing stains for CZR PRESS LF porcelain. It is necessary to use this stain after CZR PRESS LF layering. Baking temperature of this stain is 840°C.

PURE WHITE	BLUE	YELLOW	CERVICAL 1	EARTH BROWN	PINK	B+	3g each
GRAY	GREEN 1	ORANGE 1	CERVICAL 2	REDDISH BROWN	RED	C+	
BLACK	GREEN 2	ORANGE 2	CERVICAL 3	SALMON PINK	A+	D+	
GLAZE						10g	1 color
ES LIQUID						10ml	1 bottle
ES COLOR GUIDE							1 piece
TECHNICAL INSTRUCTIONS							1 copy



NORITAKE ACT LIQUID

This is Anti-Greening liquid for EX-3 Press material. Sometimes it is necessary to use ACT liquid for any alloy that contains some silver (Ag) to prevent greening.



ACT LIQUID	20ml
------------	------

EX-3 PRESS DISPO PLUNGER

NORITAKE EX-3 PRESS DISPO PLUNGER is disposable plunger set for EX-3 PRESS. 2 kinds of size are available for each 2g & 5g EX-3 PRESS INGOT.

LF LIQUID

EX-3 PRESS LF porcelain should EXCLUSIVELY be mixed with LF LIQUID.